INHALT

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1.	Intended use
2.	Scope of delivery
3.	Operating concept
4.	Performance Characteristics
5.	Qualification/Installation
6.	Setup
6.1	Preparation
6.2	Usage
7.	Adjusting the respiratory parameters
8.	Disassembly
9.	Cleaning
10.	Maintenance intervals

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For your own and you patient's safety:

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PLEASE READ THIS MANUAL VERY CAREFULLY

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Every interaction with the apparatus requires knowledge and observance of this very manual. The device is ONLY to be used according to this document.

Maintenance

The device must be inspected and maintained annually by professionals. Maintenance on this device is only to be performed by trained professionals.

For the completion of a maintenance contract, as well as for maintenance in general, we recommend our Next3D service.

ONLY original parts by Next3D are to be used for repairs. Please consult the chapter "Maintenance intervals"

Safety inspections

The device has to undergo periodical inspections. See "Maintenance intervals" for more detailed information

Accessories

Only accessories provided in our official list are to be used with the device.

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Liability regarding functionality and damages.

Non-certified medical product, only for exceptional cases, if no other certified product is available!

The user or owner of the device is in liable in case of deployment.

Liability for the function of the apparatus is transferred to the owner or operator in any case if the device is improperly maintained or repaired by persons not belonging to Next3D Service or if it is handled in a way that does not correspond to its intended use. Next3D shall not be liable for damage that occurs due to failure to observe the above instructions. Warranty and liability conditions of the Next3D terms of sale and delivery are not extended by the above notices.

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Airlog 2020 -

Volume controlled emergency ventilation system

1. INTENDED USE

Areas of deployment

The system is only intended for emergency ventilation, as part of which volume-controlled ventilation can be carried out. It does not provide assisted ventilation and was not designed to recognize spontaneous breathing and to adapt to it (no assisting device). In contrast to the available conventional and certified ventilation devices, which are differentiated with regard to the ventilation modes, it is only variable in the ventilation parameters, which are essential, but are life-supporting in the case of existing respiratory insufficiency. Self-breathing is only replaced in a volume-controlled manner.

Mobile and stationary deployment in emergency or primary care of emergency patients, while transportation in water-, air- or floor-bound rescue equipment, during relocation trips, during relocation inside of a clinic, during secondary transport from one clinic to another, in the emergency room — as long as a **permanent** power supply is present.

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Low voltage device does not work on battery power!

Observation of ventilation

While ventilating the patient the device's correct function is to supervised and checked regularly by qualified medical staff.

Keep a manual ventilation at hand

If, through a malfunction, the life preserving function of the apparatus is not ensured, the ventilation has to be guaranteed by the use of an independent ventilation device which has to be put into service immediately.

DO NOT use the device in areas endangered by explosion! Danger of explosion! DO NOT use the device in toxic areas!

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2. SCOPE OF DELIVERY

- 1x Emergency Ventilation System Airlog 2020
- 1x Sensor unit with volume and pressure sensors
- 1x Power supply unit
- 1x IFC lead

Parts **NOT** included in the package

- Ambu-Bag
- Valve assembly of the Ambu-Bag
- Peep-Valve
- Filter
- Tracheal tube
- Respiratory mask

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Display and adjustment of the following parameters

- Respiratory volume in ml (Vt)
- Respiratory frequency in breaths per minute
- Limit for inspiratory pressure for the alarm signal to sound (mbar)

Acoustic and visual alarm signal when pressure limit is exceeded or at disconnection

Ambu-Bag, tube system not included in package contents:


Consists of biocompatible, steam sterilizable synthetic material

The emergency ventilation system Airlog 2020 was tested and evaluated with the lung phantom LS800 made by Dräger. The testing system simulates the following parameters:

- Resistance 2-128 mbar/L/s
- Compliance 0,1-0,01 L/mbar.

5. QUALIFICATION / INSTALLATION

The product may only by used after first time installation and training by Next3D staff.

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6.1 Preparation

Visual inspection for damages on the system

6. SETTING UP THE AIRLOG 2020

Volume controlled ventilation for respiratory

frequencies from 4 to 25 breaths (AF) per min.

A breathing rate (AF) of 26-30 per minute and

above is not intended for the proper use of the

Airlog 2020 due to its material properties and is

solely at your own risk. Ventilation is performed

via ambient air. An additional oxygen supply is

possible via the Ambu-Bag. This must contain a

pressure dosing valve provided for this purpose.

to compressed gases is prohibited.

The direct connection of the Airlog 2020 ventilator

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
6.2 Usage

Placing the Ambu-Bag

The Ambu-Bag intended for volume-controlled ventilation and functioning properly is placed in the device provided for this purpose and the fixation bar (clamp) is closed and then locked in place with a screw (see Figure 2,3). A mechanical test confirms the stable fit / fixed position of the Ambu-Bag. The user ensures that the Ambu-Bag is regularly compressible and generates a sufficient air expulsion.

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Fixing the Ambu-Bag with the supplied clamp

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Connecting the sensor unit to the system

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Plugging in the electrical connection for the sensor assembly

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
Connecting the filtering unit to the sensor assembly

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Connecting the tracheal tube to the filter

The system must only be used as shown in the image below and with accessories shown.

Intubation of the patient. Long tubes may not be used as extension.

Connecting the power supply to the unit

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7. ADJUSTING THE RESPIRATORY PARAMETERS

Turn the system on by pressing the On/Off switch

Adjusting the respiratory volume per breath in ml

The adjustment of the RV takes places as needed and according to before measured values of $\rm O_2$ and $\rm CO_2$ partial pressure (blood/pulsoxymetric), the predominate disease, as well as to the treatment and therapeutic aims. The adjustment is to be done by trained intensive care staff and has to be carried out according to physiological and pathophysiological conditions.

Application of the ventilation gas (air mix, O2 admixture) will only occur through and in the used Ambu-Bag. The use of a ventilation gas mix via a central supply (compressed air) is prohibited. The application of bottle gas (e.g. O2) for insufflation via Ambu-Bag is only permitted by pressure reducers and only in accordance with the instructions for use of the respective manufacturer.

Adjusting the respiratory frequency

The adjustment of the RF in breaths per minute is done as needed and according to before measured values of $\rm O_2$ and $\rm CO_2$ partial pressure (blood/pulsoxymetric), the predominate disease, as well as to the treatment and therapeutic aims. The adjustment is to be done by trained intensive care staff and has to be carried out according to physiological and pathophysiological conditions.

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Setting p_{max} for the alarm function

If the patient is connected to the system:

- Respiratory volume check setting and correct value specific to patient
- 2. Check pressure value shown on LCD.
- Set upper pressure limit pmax. If the upper limit is reached the system plays an alarm signal and shows a visual indication as well. The ventilation is continued by the system. The operating staff has to act appropriately.

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8. DISASSEMBLY

Dismantling is carried out in reverse order of the steps described under point 8. The Ambu-Bag is removed after disconnecting the sensor system by loosening the fixing bracket (clamp).

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9. CLEANING

- Cleaning is done by using alcohol infested cloth(s).
- Do **not** wet clean the apparatus and protect the electric motor from moisture.
- No fluids may enter the unit.
- Only use accessories **once**, according to their specific manual.
- Dispose of single-use items according to national regulations.

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10. MAINTENANCE INTERVALS

The apparatus or its parts have to disinfected and cleaned before any kind of maintenance, even when returning the system for repairs!

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Disposal

- After retiring the Airlog 2020, consult proper enterprises and take care of correct disposal of the system.
- Obey the relevant laws of your country.

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11. PLEASE PAY ATTENTION TO THE FOLLOWING GUIDELINES!

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Manufacturer: Next 3D, Liebigstr. 13, 04103 Leipzig

The device is a non-certified medical product.

The device is not sterilized on arrival. This has to be taken care of by the owner or user.

The product may only be one-time sterilized

The apparatus is fragile and thus, is to be handled with care.

The device may only be stored at room temperature.

The device has to put into the supplied package for storage, to protect the unit from UV radiation.

Disposal may only occur under the rules of according laws.

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12. REPORTING INCIDENTS

All incidents related to or caused by the product are to be reported to the manufacturer and to the proper state authorities of the state or country the patient is living or has lived in.

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Volume controlled ventilation

Respiratory volume RV:

 50ml - 500ml (depending on the installed Ambu-Bag)

Respiratory frequency RF:

Volume controlled ventilation for respiratory frequencies from 4 to 25 breaths (AF) per min. A breathing rate (AF) of 26-30 per minute and above is not intended for the proper use of the Airlog 2020 due to its material properties and is solely at your own risk.

The emergency ventilation system Airlog 2020 was tested and evaluated with the lung phantom LS800 made by Dräger. The testing system simulates the following parameters:

- Resistance 2-128 mbar/L/s
- Compliance 0,1-0,01 L/mbar.

22 23